



GLOBUS

RF CLINIC skin

Professional Radiofrequency

11 PROGRAMS

4 Face - 7 Skin

GLOBUS RF CLINIC skin is an innovative PROFESSIONAL machine generating radiofrequency waves for aesthetic applications, against ageing of the facial skin.

Globus RF Clinic skin promote as well the Vascularisation and the regeneration of the cutaneous tissues, eliminating the skin imperfections caused by the sun exposition and the skin ageing.

RADIOFREQUENCY APPLICATIONS

Many scientific researches have demonstrated that RADIOFREQUENCY stimulate the blood circulation and the cellular metabolism by producing:

- wrinkle smoothing and skin smoothing
- recovery of cutaneous elasticity (elastosis and loose skin)

Thanks to GLOBUS Technology, it is possible to control all the impulses parameters (emission frequency from 470 KHz to 2 MHz, time and power), and to modulate the energy in the treated area in a way to obtain the maximum efficiency both in the superficial and deep treatments.

HOW DOES IT WORK?

Radiofrequency (RF) produce energy, perceived as heat, which stimulate the metabolism of the treated tissues. The radiofrequency flow, crossing the skin, increase the temperature of the lower tissues, and in particular of the supporting structure, elastin and collagen, which determine the skin elasticity and firmness.

The heat generated activates the fibroblasts pushing them to produce new collagen enhancing the tissue regeneration (collagenesis). The epidermis turns out to be more tonic and the skin tighter.

• 10 FREE PROGRAMS

• 4 programs FACE

Area dedicated to the face treatments

• 7 programs SKIN

Area dedicated to the skin from various parts of the body



TECHNICAL FEATURES

Output Power: 100 W

Output: 1 - Head: 1

Emission frequency: from 470 KHz to 2 MHz

Power Supply: main 220V 50/60 Hz and rechargeable batteries (in option)

Dimensions: 17 x 22 x 6 cm

Backlight Display

Resistive system

CE Device

Standard: EN 60601-1; EN 60601-1-2

EQUIPMENT

- 1 carrying case
- 1 unit Globus RF Clinic skin
- 1 Head RF Globus
- 1 operating manual
- 1 contact gel bottle



GLOBUS

GLOBUS ITALIA S.R.L. - Via Vittorio Veneto 52 - 31013 Codognè (TV) - ITALY

Tel (+39)0438.7933 - Fax (+39)0438.793363 - info@globuscorporation.com - www.globuscorporation.com



GLOBUS

RF CLINIC skin

Program List

FACE

Wrinkles face
Wrinkles Eyes
Acne Face
Micro-lifting face

SKIN

Wrinkles body
Wrinkles sensitive skin
Micro-Lifting décolleté
Stretch marks
Skin elasticity
Collagen Bio-Skin
Scars



GLOBUS

GLOBUS ITALIA S.R.L. - Via Vittorio Veneto 52 - 31013 Codognè (TV) - ITALY
Tel (+39)0438.7933 - Fax (+39)0438.793363 - info@globuscorporation.com - www.globuscorporation.com