

FOR EFFECTIVE PSORIASIS THERAPY



TECHNOLOGY FOR THERAPY

PSORIAMED[®] is a clinically tested device concept exactly based on the requirements of psoriasis therapy that is quickly effective and can be implemented easily through the application of stimulation current.

The special interference currents (IFC) overcome the resistance barriers of the skin without additionally irritating or damaging the affected regions, thus acting on the functionally disturbed areas. The highly antipsoriatic effect of IFC is owned to the normalization of the segmentation processes (in case of hyper proliferation) as well as to its anti-inflammatory and antioedematous effect.

Therapy with PSORIAMED[®] stimulates cellular self-help by transmission of adequate information on a cellular level with short treatments with pleasantly low current density. The highly efficient therapy, which is free of sideeffects, represents a perfect treatment alternative especially for patients resistant to therapy.

To achieve an efficient therapy, the treatment has to be performed twice a day (mornings and evenings). Application of stimulation currents is either performed under water or using special electrodes.

The device menu provides pre-programmed therapy suggestions and therapy images. The device is easy to use due to the clearly structured display and the PHYSIOMED one-button operation. For application under water a foot switch can be used to adjust the amperage without problems. Battery operation ensures safe operation of the device and enables you to use it conveniently outside the practice.

TECHNICAL DATA:

Max. power output: Power consumption: Rated voltage: Protection class: Dimensions $(W \times H \times D)$: Weight:

100 mA/500 Ohm 25 VA 100-240 V 1, type BF

315 x 104 x 285 mm 3.1 kg

STANDARD ACCESSORIES:

- 2 Treatment shells
- 2 Plate electrodes EF 400
- 2 Viscose covers EF 400
- 10 Spacer grids
- 2 Elastic velcro straps 10 x 125 cm
- 1 Patient lead (bipolar)
- 1 Mains cable
- 1 Transportation bag
- 1 Operating instructions
- 1 Therapy manual "Psoriasis"

ADDITIONAL ACCESSORIES:

Footswitch Special electrode (arm) Special electrode (head) Connection cable with clamp fastening (set of two) Plate electrode EF 200 (set of 2) Viscose cover EF 200 Trolley Universal II

PHYSIOMED devices comply with up-todate safety standards and the medical device regulations. They are manufactured according to current quality insurance standards. Our devices are tested for type examination and have the official CE marking according to regulation 93/42/EEC.

Ask for our attractive leasing offers!

e

ш

Σ

∢

2

0

S

Δ.

Clearly structured display with treatment images

TECHNOLOGY FOR THERAPY



Footswitch (optional) For safe and convenient operation during the treatment.

PHYSIOMED ELEKTROMEDIZIN AG Hutweide 10 91220 Schnaittach/Laipersdorf Germany



PHONE +49 (0) 91 26/25 87-0 FAX E-MAIL WEB

C E 1275

+49 (0) 91 26/25 87-25 info@physiomed.de www.physiomed.de