

# HiToP® 184



4 independent channels

Display LCD

Operating frequencies

Channels

Timer Supply voltage Power consumption

Protection class Equipment class Dimension

Weight Warranty 0.1 - 200 Hz 4096 – 32768 Hz 0 - 90 min115/230 V, 50/60 Hz max 100 VA I/BF (IEC 601-1) IIa (acc. to MDD)

15,5 x 42 x 40,5 cm (H x W x L)

13 kg 24 months

Subject to technical changes without prior notice.

# HiToP® 184

## A four-channel pain therapy device for a holistic treatment of the patient.

The HiToP® 184 permits treatment of various indications. The main indications are gonarthrosis, coxarthrosis, cervical spine syndrome, lumbar spine syndrome, shoulder-arm syndrome. Moreover, polyneuropathia, epicondylitis and other indications are treated successfully. The simultaneous treatment of

cause and symptom is the reason for the good results of this therapy. Parallel to the various indications the device also offers whole body vitalization.

A part from the medical field HiToP® 184 offers an extremely wide choice of programs in the wellness area. Lipolysis, cellulite, fat burning, sport, body shaping are programs that can be applied parallel to the whole body vitalization.

## **Medical programs:**


- SimulFAM® i
- Reha
- Mobilization
- Musclepump
- Trauma
- Lymphdrainage
- Edema
- Pain
- Nerveblock
- Universal

### Wellness programs:

- Vitalisation
- Sport
- Elektro massage
- Shaking massage
- Lipolysis
- Fatburning
- Bodyshaping
- Cellulite

### The most important advantages at a glance:

- 4 independent channels for the treatment in two
- Whole body treatment for vitalisation
- Easy operation by selecting programs from a list.
- Safety circuit to avoid voltage sensations caused by bad electrode placement.
- No side effects.



gbo Medizintechnik AG Kleiststrasse 6 D-64668 Rimbach

Germany

Tel.: +49(0)6253/808-0 Fax: +49(0)6253/808-245 e-mail: Info@gbo-med.de internet: www.gbo-med.de